## **Symmetry Classification for Serial Crystallography Experiments**

Groups with white backgrounds are merohedral and will exhibit indexing ambiguities. Move directly downwards to the nearest cell with a shaded background to find the corresponding "source symmetry".

Chiral groups are shown in bold, centrosymmetric groups are underlined.

		oint G	40UDC		<u> </u>	Space Croups						
		omit Gi	toups		Space Groups							
Triclinic 1	lattice		1			P4	7.1					
1 <b>1</b>					P1				P1			
Monoclin	ic lattice				_							
	2		m		P2, P2 <sub>1</sub> , C2			Pm, Pc, Cm, Cc				
			<u>2/m</u>						2/m, P2 <sub>1</sub> /m, C2/m, P2/c, P2 <sub>1</sub> /c, C2/c			
Orthorhombic lattice												
	222		mm	2	P222, P222 <sub>1</sub> , P2 <sub>1</sub> 2 <sub>1</sub> 2, P2 <sub>1</sub> 2 <sub>1</sub> 2 <sub>1</sub> , C222 <sub>1</sub> , C222, F222, I2 <sub>1</sub> 2 <sub>1</sub> 2 <sub>1</sub>			Pmm2, Pmc2 <sub>1</sub> , Pcc2, Pma2, Pca2 <sub>1</sub> , Pnc2, Pmn2 <sub>1</sub> , Pba2, Pna2 <sub>1</sub> , Pnn2, Cmm2, Cmc2 <sub>1</sub> , Ccc2, Amm2, Aem2, Ama2, Aea2, Fmm2, Fdd2, Imm2, Iba2, Ima2				
			mmm					Pmmm, Pnnn, Pccm, Pban, Pmma, Pnna, Pmna, Pcca, Pbam, Pccn, Pbcm, Pnnm, Pmmn, Pbcn, Pbca, Pnma, Cmcm, Cmce, Cmmm, Cccm, Cmme, Ccce, Fmmm, Fddd, Immm, Ibam, Ibca, Imma				
Tetragona	al lattice											
4	4 4			<u>4/m</u>	P4, P4 <sub>1</sub> , P4 <sub>2</sub> , P4 <sub>3</sub> , I4, I4 <sub>1</sub>	P4, I4				P4/m, P4 <sub>2</sub> /m, P4/n, P4 <sub>2</sub> /n, I4/m, I4 <sub>1</sub> /a		
422	42m	$0m \mid 4m0 \mid 4mm \mid \blacksquare \qquad \stackrel{\bullet}{=} \stackrel{\bullet}{=} \qquad \mid$		P42m, P42c, P42 <sub>1</sub> m, P42 <sub>1</sub> c, I42m, I42d	1	n2, P4c2, P4b2, n2, I4m2, I4c2	P4mm, P4bm, P4 <sub>2</sub> cm, P4 <sub>2</sub> nm, P4cc, P4nc, P4 <sub>2</sub> mc, P4 <sub>2</sub> bc, I4mm, I4cm, I4 <sub>1</sub> md, I4 <sub>1</sub> cd					
<u>4/mmm</u>						P4/mmm, P4/mcc, P4/nbm, P4/nnc, P4/mbm, P4/mnc, P4/nmm, P4/ncc, P4 <sub>2</sub> /mcm, P4 <sub>2</sub> /nbc, P4 <sub>2</sub> /nbc, P4 <sub>2</sub> /mnm, P4 <sub>2</sub> /nmc, P4 <sub>2</sub> /ncm, I4/mmm, I4/mcm, I4 <sub>1</sub> /amd, I4 <sub>1</sub> /acd						

Trigonal lattice

3	<u>3</u>		R3 (H3)	<u>R3 (H3)</u>
32		3m	R32 (H32)	R3m (H3m), R3c (H3c)
		<u>8m</u>		$\overline{R3m}$ ( $\overline{H3m}$ ), $\overline{R3c}$ ( $\overline{H3c}$ )

Hexagonal lattice

1102	Hexagonal lattice																
	3	3		<u>1</u>	<u>3</u>					<b>P3</b> , <b>P</b> 3	3 <sub>1</sub> , P3 <sub>2</sub>			<u>P3</u>			
			3:	m1	31	lm						P3m1	l, P3c1	P31m	ı, P31c		
									P6, P6 <sub>1</sub> ,				,		,		
6					<u>5</u>			<u>6/m</u>	P6 <sub>5</sub> , P6 <sub>2</sub> ,					P <del></del> 6			<u>P6/m</u> ,
									P6 <sub>4</sub> , P6 <sub>3</sub>								<u>P6<sub>3</sub>/m</u>
	312					<u>31m</u>				P312, P3 <sub>1</sub> 12, P3 <sub>2</sub> 12					P31m, P31c		
		321	<u>3</u> m1								P321, P3 <sub>1</sub> 21, P3 <sub>2</sub> 21	<u>P3m1</u> , <u>P3c1</u>	_				
				6m2	<del>6</del> 2m		6mm						P6m2, P6c2	2 P <del>6</del> 2m, P <del>6</del> 2c		P6mm, P6cc, P6 <sub>3</sub> cm, P6 <sub>3</sub> mc	
	622	622 <u>6/mmm</u>					P622, P6 <sub>1</sub>	P6 <sub>4</sub> 22, P6 <sub>5</sub> 22, P6 <sub>2</sub> 22, P6 <sub>4</sub> 22, P6 <sub>3</sub> 22  P6/mmm, P6/mcc, P6			<u>P6<sub>3</sub>/mcm</u> , <u>P6</u>	<u>3/mmc</u>					

## Cubic lattice

23		<u>m</u> 3	P23, F23, I23, P2 <sub>1</sub> 3, I2 <sub>1</sub> 3		<u>Pm3</u> , <u>Pn3</u> , <u>Fm3</u> , <u>Fd3</u> , <u>Im3</u> , <u>Pa3</u> , <u>Ia3</u>
432	<del>4</del> 32		P432, P4 <sub>2</sub> 32, F432, F4 <sub>1</sub> 32, I432, P4 <sub>3</sub> 32, P4 <sub>1</sub> 32, I4 <sub>1</sub> 32	P43m, F43m, I43m, P43n, F43c, I43d	
	m <del>3</del> m				<u>Pm3m</u> , <u>Pn3n</u> , <u>Pm3n</u> , <u>Pn3m</u> , <u>Fm3m</u> , <u>Fm3c</u> , <u>Fd3m</u> , <u>Fd3c</u> , <u>Im3m</u> , <u>Ia3d</u>